

Care of Unit:

The Gooseneck & Super Clamp can be wiped clean with any household multi-purpose cleaner and disinfectant. We recommend Simple Green, which is a non toxic biodegradable all-purpose cleaner.

Do not submerge the unit, as it will damage the electrical components.

Do not use abrasive cleaners, as they will scratch the surface of the unit.



For Technical Support:

Call our Technical Service Department
Monday through Friday, 9 a.m. to 5 p.m. (EST)
1-800-832-8697

customer_support@enablingdevices.com

Gooseneck & Super Clamp #1549

USER'S GUIDE



50 Broadway
Hawthorne, NY 10532
Tel. 914.747.3070 / Fax 914.747.3480
Toll Free 800.832.8697
www.enablingdevices.com

Allows Convenient Positioning

Now you can position a switch anywhere you'd like with our rugged mounting system. Equipped with a flexible, extra stiff gooseneck, you can attach it to a table or wheelchair using a unique clamp. Your switch can be attached to the gooseneck with any standard mounting plate. *Size: 20¹/₄"L. Weight: 2¹/₄ lbs.*

Operation:

1. **Clamping Options:** Our mounting clamp is designed for mounting to tubular (wheelchair, bedrail, walker, etc.) or flat surfaces such as a table top or desk. For flat surface mounting use the included wedge insert shoe stored in front of the locking knob on the clamp (see Fig.1 & 2) Loosen locking knob slightly to remove wedge, tighten knob once wedge has been removed.



Fig.1



Fig. 2

2. Using the mounting clamp, attach the Gooseneck to a Wheel Chair, Table or Desk. Tighten the clamp enough to hold the Gooseneck securely, but do not over-Tighten the clamp.



Fig.3

3. To connect your mounting plate to the Gooseneck, your mounting plate must have a standard 1/4-20 threaded post to screw into the top of the Gooseneck. Our #7808 Mounting Plate for Jumbo Switch, #7809 Mounting Plate for Gumball Switch, #2080 Circular, #2082 Rectangular Mounting plates all use the standard post and are compatible with this Gooseneck Mount.
4. To mount screw your plate on until flush and tight with the top of the Gooseneck. See Fig.4.



Fig. 4