
Lighted Vibrating Mirror #358-M

USER'S GUIDE



For Technical Support:
Call our Technical Service Department
Monday through Friday, 9 a.m. to 5 p.m. (EST)
1-800-832-8697
customer_support@enablingdevices.com



50 Broadway
Hawthorne, NY 10532
Tel. 914.747.3070 / Fax 914.747.3480
Toll Free 800.832.8697
www.enablingdevices.com

Lights, vibration, action!

Use our vibrating mirror during circle time activities. This mirror encourages special explorations. When you pick up the mirror the lights turn on while the unit vibrates. This encourages the user to visually attend to himself and increases bi-lateral hand usage. The level of vibration is adjustable. Size: 8½"L x 12"W. Requires 6 AA Batteries. Weight: 1¾ lbs.

Operation:

1. The Lighted Vibrating Mirror requires 6 AA batteries. Alkaline batteries are recommended because they last longer. Do not use rechargeable batteries because they supply lower voltage and the unit may not perform as well.
2. The battery compartments are located under the unit. Carefully turn the unit over. The battery case covers are secured with one small screw each. Remove the screws by turning them counterclockwise with a Phillips head screwdriver. Install 2 AA batteries in each of the compartments, being careful to observe proper battery polarity.
3. Replace the battery case covers and secure with the screws. Do not over-tighten.
4. The ON/OFF/VIBRATION control knob is located next to the battery compartments. Turning this knob clockwise turns the unit on as well as control the vibration intensity. Turning clockwise increases vibration, counterclockwise decreases and turns the unit off.
5. To activate the Lighted Vibrating Mirror lift the unit from a lying or horizontal position and bring it to a vertical position. This will reward the user with lights and vibration. Putting the Lighted Vibrating Mirror down in a lying position will stop the lights and vibration.
6. When the Lighted Vibrating Mirror is not in use please make sure that the ON/OFF/VIBRATION control knob is in the off position. This will maximize battery life. When storing the Lighted Vibrating Mirror for long periods of time we recommend removing the batteries.

Troubleshooting:

Problem: There are no light or vibration.

Action #1: Rotate (spin) the batteries within the battery compartment to ensure a good connection.

Action #2: Replace batteries paying close attention to proper battery polarity.

Action #3: Check all batteries for freshness, if weak or dead replace with new batteries.

Action #4: Make sure that the bead chains are making contact with the metal contact rod described above, when they are moving.

Care of Unit:

The Lighted Vibrating Mirror can be wiped clean with any household multi-purpose cleaner and disinfectant.

Do not submerge the unit, as this could damage the electrical components.

Do not use abrasive cleaners, as they will scratch the surface of the unit.