# Pull Ball #416

**USER'S GUIDE** 



For Technical Support: Call our Technical Service Department Monday through Friday, 9 a.m. to 5 p.m. (EST) 1-800-832-8697 customer\_support@enablingdevices.com



50 Broadway Hawthorne, NY 10532 Tel. 914.747.3070 / Fax 914.747.3480 Toll Free 800.832.8697 www.enablingdevices.com

# Easy fun with melody and illumination!

Gentle tugging on this multifaceted spring ball will activate music and lights. The tetherball will fly back into its case when released. Increases eye-hand coordination, fine motor skills and finger isolation. To avoid cord rupture, gentle tugging of ball is recommended. Size:  $12\frac{1}{2}$ "L x 4"W x 7<sup>1</sup>/<sub>2</sub>"H. Requires 2 AA Batteries. Weight: 1<sup>1</sup>/<sub>4</sub> lbs.

#### **Operation:**

- 1. The Pull Ball requires 2 AA batteries to operate. Alkaline batteries are recommended because they last longer. Do not use rechargeable batteries because they supply lower voltage and the unit may not perform as well.
- 2. The battery compartment is located on the bottom of the unit. The compartment's cover is secured with a small screw as a safety precaution. Use a Philips head screwdriver to remove the small screws and slide the battery compartment cover off. Install the 2 AA batteries paying careful attention to proper (+) and (-) battery polarity.
- 3. Replace the battery compartment cover and re-secure the screw.
- 4. When you pull the ball in any way, you will be rewarded with music and lights!
- 5. Once you stop pulling the ball, there will be a brief delay, and then the music and lights will stop.

## **Troubleshooting:**

**Problem:** The unit does not operate, or operates incorrectly

**Action:** Check all batteries for freshness and proper (+) and (-) polarity.

**Problem:** One or more of the features is not operating, while other features are operating.

**Action:** Check <u>ALL</u> of the batteries for <u>freshness</u> and <u>proper polarity</u>. Replace any weak/dead batteries paying close attention to proper polarity.

### Care of Unit:

The Pull Ball can be wiped clean with any household multi-purpose cleaner and disinfectant.

**Do not submerge** the unit, as it will damage the electrical components.

**Do not use abrasive cleaners**, as they will scratch the surface of the unit.