

---

---

## Rock N Roll Stuart #8083

---

---

### USER'S GUIDE



**For Technical Support:**

Call our Technical Service Department  
Monday through Friday, 9 a.m. to 5 p.m. (EST)  
1-800-832-8697

[customer\\_support@enablingdevices.com](mailto:customer_support@enablingdevices.com)



50 Broadway  
Hawthorne, NY 10532  
Tel. 914.747.3070 / Fax 914.747.3480  
Toll Free 800.832.8697  
[www.enablingdevices.com](http://www.enablingdevices.com)

### **From Despicable Me!**

Activate your switch and Minion Stuart “plays” his electric guitar and moves to the beat of the music. We’ve included a sturdy stand with suction cups for easy play. *Size: 6”L x 5”W x 9½”H. Requires 4 AA Batteries. Weight: 2 lbs.*

#### **Operation:**

1. Stuart require 4 AA batteries (Not included) to operate. Use alkaline batteries only (e.g. Duracell or Energizer brand). Do not use rechargeable batteries or any other type of batteries because they supply lower voltage and the unit will not perform properly. Never mix old and new batteries together or different brands or types together.
2. The battery compartment is located is the back of the doll. Locate the Velcro opening and then using a Phillips head screwdriver open and remove the batter cover. Next install 4 AA batteries paying close attention to proper battery polarity.
3. Replace the battery compartment cover and secure with screw. Set the TRY ME OFF PLAY switch to the PLAY position.
4. Connect Stuart to any external capability switch (not included) via the jack located at the end of their cord. Stuart is equipped with a 1/8" jack for use with external capability switches with a 1/8" plug. Make sure to plug your switch all the way in. There should be no gaps.
5. Actuate your external capability switch or press the top of his left foot, he will “play” his electric guitar and moves to the beat of the music. Each activation of your switch or Stuart’s paw will make him say something new.

6. We have also included a sturdy stand to hold Stuart upright for easy play. Slide Stuart’s Guitar strap on the back over the upright part of the stand to help keep him in place while using the stand.
7. When finished using Stuart turn his ON/OFF switch to off to prevent false activations and to prevent battery drain.

#### **Troubleshooting:**

**Problem:** Stuart does not function correctly.

**Action #1:** Replace the batteries paying close attention to proper battery polarity.

**Action #2:** Make sure that the connection between the external capability switch and the jack is pushed in all the way. There should be NO gaps.

**Action #3:** Try your switch with another switch operated toy/device in an effort to rule out your switch as the source of the problems.

#### **Care of Unit:**

If you will not be using the Stuart for an extended amount of time, we recommend removing the batteries.

**Do not submerge** the unit in water, as it will damage the electrical components.

Stuart is surface washable only.