
**Vibrating Animals #9405 White Rabbit,
Basset Hound #9406 #9319 Floppy Rabbit**

USER'S GUIDE



For Technical Support:

Call our Technical Service Department
Monday through Friday, 9 a.m. to 5 p.m. (EST)
1-800-832-8697

customer_support@enablingdevices.com



50 Broadway
Hawthorne, NY 10532
Tel. 914.747.3070 / Fax 914.747.3480
Toll Free 800.832.8697
www.enablingdevices.com

Great tactile awareness!

Give a hug and these soft, cuddly animals will gently vibrate, providing somatosensory stimulation and relaxation. These vibrating animals will be friends for life. Size: 8" to 20" long depending on selection. Requires 2 AA Batteries. Weight: ¾ lb.

Operation:

1. Vibrating Animal requires 2 AA Batteries to operate. Alkaline batteries must be used (e.g. Duracell or Energizer Brands). Do not use rechargeable batteries or any other type of batteries because they supply lower voltage and the unit will not perform properly. Never mix old and new batteries together or different brands or types together.
2. The battery compartment is located inside Vibrating animal's lower pouch opening. Remove the Velcro attached vibration unit. Next slide the battery cover off and install 2 AA batteries paying close attention to proper battery polarity. Replace the battery compartment cover.
3. Replace the vibration unit back inside the Animal making sure the Velcro side of the vibration unit secures to the Velcro strip inside the Animal. Each time you squeeze or hug the Animal you will feel vibration.

Troubleshooting:

Problem: Vibrating Animal does not seem to function properly.

Action: Low battery power may affect toy function, including irregular or intermittent vibration. Replace all batteries with new alkaline batteries.

Care of Unit:

If you will not be using Vibrating Animal for an extended amount of time, we recommend removing the batteries.

Do not submerge the unit in water, as it will damage the electrical components. Vibrating Animal are surface washable with cold water, air dry only.