Vibrating Forest Friends #3934 (Set of 4) Baby Rabbit #3930, Bear #3932 & Turtle #3933

USER'S GUIDE



For Technical Support: Call our Technical Service Department Monday through Friday, 9 a.m. to 5 p.m. (EST) 1-800-832-8697

customer_support@enablingdevices.com



50 Broadway Hawthorne, NY 10532 Tel. 914.747.3070 / Fax 914.747.3480 Toll Free 800.832.8697 www.enablingdevices.com

Soothing vibration!

Give one of these soft plush Forest Friends a hug and feel gentle vibration. They are fantastic for providing proprioceptive input, tactile stimulation, and cause and effect. Each requires 2 AA Batteries. Weight: $\frac{1}{2}$ lb. Size: 10"H x 5"W x 5"D.

Operation:

- Vibrating Forest Friends requires 2 AA Batteries to operate. Alkaline batteries must be used (e.g. Duracell or Energizer Brands. Do not use rechargeable batteries or any other type of batteries because they supply lower voltage and the unit will not perform properly. Never mix old and new batteries together or different brands or types together.
- 2. The battery compartment is located inside Vibrating Forest Friend's lower pouch opening. Remove the Velcro attached vibration unit. Next slide the battery cover off and install 2 AA batteries paying close attention to proper battery polarity. Replace the battery compartment cover.
- 3. Replace the vibration unit back inside the Forest Friend making sure the Velcro side of the vibration unit secures to the Velcro strip inside the Forest Friend. Each time you squeeze or hug the Forest Friend you will feel vibration.

Troubleshooting:

Problem: Vibrating Forest Friends does not seem to function properly.

Action: Low battery power may affect toy function, including irregular or intermittent vibration. Replace all batteries with new alkaline batteries.

Care of Unit:

If you will not be using Vibrating Forest Friends for an extended amount of time, we recommend removing the batteries.

Do not submerge the unit in water, as it will damage the electrical components Vibrating Forest Friends are surface washable with cold water, air dry only.