
Vibrating Sensory Pod #1169

USER'S GUIDE



For Technical Support:

Call our Technical Service Department
Monday through Friday, 9 a.m. to 5 p.m. (EST)
1-800-832-8697

customer_support@enablingdevices.com



50 Broadway
Hawthorne, NY 10532
Tel. 914.747.3070 / Fax 914.747.3480
Toll Free 800.832.8697
www.enablingdevices.com

Soft sensory feedback!

Soft pillow that vibrates in your hand when you turn your wrist. Great for teaching cause and effect, developing grasping and supination skills, and providing soothing somatosensory feedback. Size: 3"D x 2"H. Requires 1 N Battery Not Included. Weight: 5 oz.

Operation:

1. The Vibrating Sensory Pod requires one N-cell battery. To install the battery, unzip the pillow pouch and remove the vibration module. Next install the N-cell battery following proper (+) and (-) polarity. Replace vibrating module and zip pouch closed.
2. The Soft pillow will vibrates in your hand when you turn your wrist. Great for teaching cause and effect, developing grasping and supination skills, and providing soothing somatosensory feedback.
3. When not in use or being stored for a prolong period of time, e.g. one month or longer, we recommend removing the batteries from the Vibrating Sensory Pod during storage to prevent battery drain or leakage, which can cause damage to the Vibrating Sensory Pod.

Troubleshooting:

Problem: The Vibrating Sensory Pod does not vibrate.

Action #1: Check that the battery is installed properly following proper (+) and (-) polarity.

Action #2: Check battery for strength, if weak or dead replace with new battery.

Care of Unit:

The Vibrating Sensory Pod pillow pouch is surface washable only. Air dry only.

Do not submerge the unit, as it will damage the contents and the electrical components.

Do not use abrasive cleaners, as they will scratch the surface of the unit.