

Bluetooth Sensory Tube - #2292



enabling
devices

50 Broadway

Hawthorne, NY 10532

Tel. 914.747.3070 / Fax 914.747.3480

Toll Free 800.832.8697

www.enablingdevices.com

Plays your music!

Red, blue, green, and multi-colored lights glow through this unique frosted tube. Each time you push a colored button, the light changes to match the button. A Bluetooth receiver allows you to play music from your iPad, iPhone, iPod Touch, Android or any Bluetooth device. Or you can skip the Bluetooth and plug your music device directly into the Mood Tube using a stereo cord (included). You choose the desired atmosphere with your music – create a soothing, relaxing ambiance or a stimulating upbeat experience. Features:

- Great for Tabletop Use
- Wireless Bluetooth
- 4 different time modes for lights – 15 seconds, 1 minute, 3 minutes or continuous
- Uses 9 volts of power

Size: 27"H x 16"Diameter. Weight: 5 lbs.

Setup:

You will first need to place your Bluetooth Sensory Tube in the location where it will be used; these factors must be present for proper operation. One the Bluetooth Sensory Tube must be placed on a flat secure level surface, e.g. floor or very sturdy table, do not place the Bluetooth Sensory Tube on any type of folding table or unstable, unlevel surface. Second the Bluetooth Sensory Tube requires a standard 120 Volt wall outlet AC power source. **Please Note:** If placing the Bluetooth Sensory Tube on the floor the power cord and any extension cords if used should be taped down to prevent any of it's users from tripping on the cords. **Please Note:** Make sure the ON/OFF switch on the Bluetooth Sensory Tube is set to OFF before plugging it into an AC wall outlet.

Setting Up Button Color Time Modes:

1. 4 different time modes
2. Push blue and Green for 15 Seconds (Default on Power Up)
3. Push Blue and red for 1 Minute timer and chose a color.
4. Push blue and Multi-color for 3 Minute timer and chose a color.
5. Push green and red to disable timer and play continuously.

Note: When in continuous mode just turn the unit off and then back on to re-enable the 15 second timer.

When turning the unit on, it defaults to 15 seconds on each color. When pressing buttons, the tube should flash to indicate change of timing. See Photo #1.



Photo #1

Bluetooth Syncing your Audio Device:

To setup your Bluetooth enabled device with the Bluetooth Sensory Tube, first turn the speaker on. The speaker has its's own power button independent of the Sensory Tube's main power button. The speaker must be turned ON/OFF separately for each use. Next follow the included Bluetooth syncing instructions. Or view online here:

https://d2c6jjk3vnoatm.cloudfront.net/spree/products/accessoies/215/8552C73DA373CA1_A7910_manual.pdf

Standard Audio Connection:

Included is an 8ft. 1/8" stereo male to male audio cable. This will allow you to connect any audio source with an earphone jack or line out jack directly to the Mood Tube and use its speakers to play back your music. This includes CD players, Radios and Tape players. You will need to disconnect the audio input jack from the Bluetooth Sensory Tube's Bluetooth Receiver and plug one end of the included audio cable into that jack and the other end into the ear phone or line out jack on your audio device. **See photos 2, 3 & 4.**



Photo #2



Photo #3



Photo #4

Troubleshooting:

If the unit is not functioning as it is supposed to:

- Check the AC power cord is plugged in the outlet all the way. Try unplugging the cord and then re-plug it back in.
- Make sure the Bluetooth speaker is turned ON.
- Check your Bluetooth setting in your device/try re-syncing to the Bluetooth Sensory Tube.
- Make sure that the AC wall outlet is working.
- Check all your settings.
- Make sure all the colored switches are moving freely.
- Review our setup video for Bluetooth pairing at:
http://enablingdevices.com/catalog/Sensory_Space/Bubble_Tubes/mood-tube
- If you have tried ALL of the above contact our Tech Support Department at customer_support@enablingdevices.com or 1-800-832-8697 Ext. 320.



For Technical Support:

Call our Technical Service Department
Monday through Friday, 9 a.m. to 5 p.m. (EST)
1-800-832-8697
customer_support@enablingdevices.com

Rev: 8-1-14

