

50 Broadway Hawthorne, NY 10532 Tel. 914.747.3070 / Fax 914.747.3480 Toll Free 800.832.8697 www.enablingdevices.com



For Technical Support: Call our Technical Service Department Monday through Friday, 9 a.m. to 5 p.m. (EST) 1-800-832-8697 customer support@enablingdevices.com

Tactile feedback!

This textured communicator provides the user with a fantastic tactile experience when they initiate language. It is great for quick communication. The different colored switches provide the user with good visual cues for choice making.

- 20 seconds total recording time
- Records 2 messages
- Built-in icon holders
- Two output jacks

Size: $8\frac{1}{4}$ "L x $2\frac{1}{4}$ "H x 7"W. Requires 4 AAA Batteries. Weight: 1 lb.

Operation:

- Carefully turn unit over and remove the small screw from the battery compartment. Slide the battery compartment cover off. Observing polarity, install four AAA size batteries (Not included) into the holder. Replace cover and screw. Alkaline batteries must be used only (e.g. Duracell or Energizer brand). Do not use rechargeable batteries or any other type of batteries because they supply lower voltage and the unit will not perform properly.
- 2. Never mix brands, types or old and new batteries together. Replace the battery compartment cover and replace the screw.
- 3. Turn the unit on by turning the ON/OFF/VOLUME knob located on the bottom panel of the unit. To record, first press and hold the red "RECORD" button press down one of the colored textured switch plates. While continuing to hold the plate down, speak your message. Release both when you finish recording. Record on the remaining switch plate in the same manner.
- 4. To play back your message, press lightly on the colored textured switch plate of your choice. It is not necessary to continue holding the switch plate down when listening to the message. The unit will continue to playback the entire message, even if pressure on the switch plate is released.
- 5. The unit will retain recorded information indefinitely, even if the batteries are removed. Previously stored information will be erased only when a new recording is performed.
- 6. Playback volume can be adjusted using the "ON/OFF/VOLUME" control knob located on the bottom panel of the unit.

To use as a switch:

- 7. This unit can also be used to activate switch adapted devices. Simply attach the device to the corresponding 1/8 inch female jack using the cords provided. There is an output jack for each or the switch plates. When the plate is pressed, the message will play and the device will be activated. Make sure all connections are plugged in all the way with no gaps between plug and device.
- 8. The Talkable has holders for placing lcons or Pictures. On the back of the unit there is a holding strip where lcons or Photos can be placed. When inserting icon in the back holding strip, insert your icon by starting at one side of the strip. Put the icon in at a slight angle for easy installation.

Please Note: When using the holding strip if you're lcons, Photos or words are printed on regular printer paper, we recommend laminating them. If not they will not stand upright or be secure in the holding strip.

Troubleshooting:

Problem: The unit fails to operate, or operates incorrectly.

Action **#1:** Make sure there is no dirt or obstacles between the plates and the housing of the unit.

Action #2: Disconnect the power by removing one of the batteries for 20 seconds in order to reset the unit. Re-install the batteries. If the unit does not work correctly, replace all the batteries.

Problem: My device connected to the Talkable does not work.

Action **#1:** Make sure there are no gaps between your cord connections. Action#2: Try your device with another switch you know works to rule it out as the source of the problem.

Care of the Unit:

The Talkable can be wiped clean with any mild household multi-purpose, nonabrasive cleaner and disinfectant. We recommend Simple Green, which is a non toxic biodegradable all-purpose cleaner.

Do not use abrasive cleaners, as they will scratch the surface of the unit.

Do not submerge the unit, as it will damage the electrical components.