

---

---

## **Somatosensory Bead Chain with Mirror #356**

---

---

### **USER'S GUIDE**



**For Technical Support:**

Call our Technical Service Department  
Monday through Friday, 9 a.m. to 5 p.m. (EST)

1-800-832-8697

[customer\\_support@enablingdevices.com](mailto:customer_support@enablingdevices.com)



50 Broadway

Hawthorne, NY 10532

Tel. 914.747.3070 / Fax 914.747.3480

Toll Free 800.832.8697

[www.enablingdevices.com](http://www.enablingdevices.com)

**Comes with a peek-a-boo mirror!**

Push aside the bead chain, and you'll see yourself in the unbreakable plastic mirror. Great for providing visual stimulation and body identification. Size: 17¼"L x 8¼"W x 13"H. No Batteries required. Weight: 3¾ lbs. Appropriate for individuals of any age.

**Care of Unit:**

The Somatosensory Musical Bead Chain/Switch can be wiped clean with any household multi-purpose cleaner and disinfectant, we recommend Simple Green, which is a non toxic biodegradable all-purpose cleaner. If the bead chain becomes dull over time, it can be wiped clean with WD-40. Care must be taken, however, to remove all residue of the oil with a multi-purpose cleaner or disinfectant.

**Do not submerge** the unit, as this could damage the electrical components.

**Do not use abrasive cleaners**, as they will scratch the surface of the unit.