

---

---

**Lights, Vibration & Music  
Somatosensory Tube #402**

---

---

**USER'S GUIDE**



**For Technical Support:**

Call our Technical Service Department  
Monday through Friday, 9 a.m. to 5 p.m. (EST)  
1-800-832-8697  
[customer\\_support@enablingdevices.com](mailto:customer_support@enablingdevices.com)



50 Broadway  
Hawthorne, NY 10532  
Tel. 914.747.3070 / Fax 914.747.3480  
Toll Free 800.832.8697  
[www.enablingdevices.com](http://www.enablingdevices.com)

### Instant rewards!

Perfect for developing grasping skills, localizing sound and teaching cause and effect. Simply pick up the tube (#3035 and #3035), hold it in a vertical position and it instantly rewards you with glowing lights and soothing vibrations. Our version (#402) adds music. All use reliable solid-state technology. Size: 2½"D x 12"L. Requires 2 AA Batteries. Weight: ¾ lb.

### OPERATING INSTRUCTIONS:

1. Remove the vinyl end cap with the small hole on the top, by holding the tip of the cap and gently pulling and turning. Only one end cap will come off, exposing the battery compartment and speed controller.
2. Pull the battery sleeve out of the tube by round cutout in its base.
3. Remove the black battery case cover by unscrewing the Phillips head screw, then slide the cover off and install two AA batteries following proper (+) and (-) markings. Replace cove and screw and slide sleeve back in and reinstall the vinyl cap.
4. While holding the tube at an upright angle, adjust the vibration speed control switch next to the round cutout starting from right to left, settings are as follows: HIGH/LOW/OFF vibration.
5. Reinstall the vinyl cover over the end of the tube.  
**Please Note: when re-installing the battery sleeve make sure to slide it just until it lines up with the end of the clear tubing.**
6. Simply turn the **Lights, Vibration & Music Sensory Tube** to a vertical position to enjoy the soothing vibrations and glowing lights.

### IMPORTANT NOTES:

If the unit fails to operate, or operates incorrectly, the following steps should be performed:

- Check batteries to ensure proper installation. If the unit still does not work correctly after the batteries are reconnected, replace all the batteries.
- **The lights on the tube do not shut off immediately, rather they fade out slowly once the tube is placed horizontally.**

If you have any problems or questions about this unit or any of our products, please call our Technical Assistance Department at 800-832-8697.