
Musical Bead Chain with Lights and Vibration

A sensory sensation!

Everyone will respond to our popular bead chain curtains. We believe you will find even greater benefits with our multi-sensory mirrored attention getters. Both versions have music and a mirror. #659 also has lights and vibration. Perfect for teaching cause and effect. Also provides fantastic auditory awareness and encourages body identification skills. #657 can also be used as a switch to activate toys or devices. Size: 17¼"L x 8¼"W x 13"H. Weight: 4¾ lbs.

Operation:

1. The battery compartments are located under the unit. Remove the screws on the battery compartment covers and slide the covers off. Observing polarity, install two AA batteries into each compartment (Not included). Spin the batteries slightly to ensure a good connection. Alkaline batteries are recommended because they last longer. Do not use rechargeable batteries because they supply lower voltage and the unit may not operate as well. Replace the battery compartment covers.

Care of Unit:

The Musical Bead Chain with Lights and Vibration can be wiped clean with any household multi-purpose cleaner and disinfectant.

If the bead chain becomes dull over time, it can be wiped clean with WD-40. Care must be taken, however, to remove all residue of the oil with a multi-purpose cleaner or disinfectant.

Do not submerge the unit.

Do not use abrasive cleaners, as they will scratch the surface of the unit.

**Musical Bead Chain w/Mirror, Light
& Vibration #659**

USER'S GUIDE



For Technical Support:

Call our Technical Service Department
Monday through Friday, 9 a.m. to 5 p.m. (EST)

1-800-832-8697

customer_support@enablingdevices.com



50 Broadway
Hawthorne, NY 10532
Tel. 914.747.3070 / Fax 914.747.3480

Toll Free 800.832.8697

www.enablingdevices.com