
Vibrating Rabbit # 9319

USER'S GUIDE



For Technical Support:

Call our Technical Service Department
Monday through Friday, 9 a.m. to 5 p.m. (EST)
1-800-832-8697

customer_support@enablingdevices.com



50 Broadway
Hawthorne, NY 10532
Tel. 914.747.3070 / Fax 914.747.3480
Toll Free 800.832.8697
www.enablingdevices.com

Great tactile awareness!

Hug this adorable seal and feel gentle vibration. Great for providing somatosensory stimulation and relaxation. *Size: 20"L. Requires 2 AA Batteries. Weight: ¾ lb.*

Operation:

1. Vibrating Rabbit requires 2 AA Batteries to operate. Alkaline batteries must be used (e.g. Duracell or Energizer Brands. Do not use rechargeable batteries or any other type of batteries because they supply lower voltage and the unit will not perform properly. Never mix old and new batteries together or different brands or types together.
2. The battery compartment is located inside Vibrating Rabbit's lower pouch opening. Remove the Velcro attached vibration unit. Next slide the battery cover off and install 2 AA batteries paying close attention to proper battery polarity. Replace the battery compartment cover.
3. Replace the vibration unit back inside the Rabbit making sure the Velcro side of the vibration unit secures to the Velcro strip inside the Rabbit. Each time you squeeze or hug the Rabbit you will feel vibration.

Troubleshooting:

Problem: Vibrating Rabbit does not seem to function properly.

Action: Low battery power may affect toy function, including irregular or intermittent vibration. Replace all batteries with new alkaline batteries.

Care of Unit:

If you will not be using Vibrating Rabbit for an extended amount of time, we recommend removing the batteries.

Do not submerge the unit in water, as it will damage the electrical components. Vibrating Rabbit are surface washable with cold water, air dry only.