
Super Sounds Guitar #1664

USER'S GUIDE



For Technical Support:

Call our Technical Service Department
Monday through Friday, 9 a.m. to 5 p.m. (EST)
1-800-832-8697

customer_support@enablingdevices.com



50 Broadway
Hawthorne, NY 10532
Tel. 914.747.3070 / Fax 914.747.3480
Toll Free 800.832.8697
www.enablingdevices.com

Open up the world of music with the Super Sounds Guitar. Adapted for maximum accessibility, the guitar adds two accessible switches and a base with suction cups for stability on a table or wheelchair tray. Practice fine motor skills by exploring the guitar's touch-sensitive neck—lightly press the cheerful animal buttons on the neck for single notes, or strum the body for a full chord. One adaptive switch activates the song mode, changing songs with each press, and the other activates the whammy bar's sound effects. Two play modes double the fun with music and silly sounds. Size: 16½"L × 8"H × 7"D. Requires 3 AA Batteries. Weight: 2 lbs.

Operation:

1. Super Sounds Guitar requires 3 AA batteries. Use alkaline batteries only (e.g. Duracell or Energizer brand). Do not use rechargeable batteries or any other type of batteries because they supply lower voltage and the unit will not perform properly. Never mix old and new batteries together or different brands or types of batteries together.
2. The battery compartment is located on the back of the Guitar. Carefully remove the small screw from the cover using a Phillips head screwdriver, and lift the cover off. Observing proper (+) & (-) battery polarity, install 3 AA size batteries into the holder. Replace the cover and screw. Next set the volume switch to low or high. Set to one of two play modes, one music or two animal sounds via the slide switch on the Guitar.
3. Activate any one of the large gumball switches on the front of the base to play fun music and sound effects. We have switch adapted two functions on the Super Sounds Guitar. The blue Einstein button at the top of the Guitar and the Whammy Bar. Our large easy to press switches will replicate these two functions as noted.
4. For a more detailed overview and information on all the functions and features, please reference the included OEM instruction guide.

Troubleshooting:

Problem: Super Star Guitar does not work.

Action #1: Rotate (spin) the batteries within the battery compartment to ensure a good connection. Check that all batteries are installed properly.

Action #2: Check all batteries for freshness, if weak or dead replace with new batteries. Use an alkaline battery (E.G. Duracell or Energizer brand).

Care of the Unit:

Super Sounds Guitar can be wiped clean with any mild household multi-purpose, non-abrasive cleaner and disinfectant. We recommend Simple Green, which is a non toxic biodegradable all-purpose cleaner.

Do not use abrasive cleaners, as they will scratch the surface of the unit.

Do not submerge the unit, as it will damage the electrical components.