Mini UV Bead Chain # 3109 (SM) & 3109R (Rock Music)

USER'S GUIDE



For Technical Support:

Call our Technical Service Department Monday through Friday, 9 a.m. to 5 p.m. (EST) 1-800-832-8697

customer support@enablingdevices.com



50 Broadway
Hawthorne, NY 10532
Tel. 914.747.3070 / Fax 914.747.3480
Toll Free 800.832.8697
www.enablingdevices.com

Sized for your tray or desktop!

Press the mirrored background and ultraviolet LEDs create dramatic colors on the colorful UV reactive bead strands. Great visual and tactile effects. Comes in two versions: one with rock music, the other with soothing music. Size: $9\frac{1}{2}$ "W x $4\frac{1}{2}$ "H x 7"D. Requires 4 AA Batteries. Weight: $1\frac{1}{4}$ lbs.

- 1. The Mini UV Mirror Bead Chain requires 4 AA batteries (not included). Alkaline batteries must only be used (e.g. Duracell or Energizer brand). Do not use rechargeable batteries or any other type of batteries because they supply lower voltage and the unit will not perform properly. The battery compartments are located under the unit's base. Gently turn the unit over to reveal the battery compartments. Remove the small screw from each cover using a Phillips head screwdriver, and then lift the battery compartment covers off. Install batteries into the holder. When installing new batteries, always be careful to observe battery polarity. Replace battery covers and screws. Never mix old and new batteries together or different brands or types together.
- 2. Press the mirrored background and ultraviolet LEDs create dramatic colors on the colorful UV reactive bead strands. Great visual and tactile effects. Release pressure from mirrored background and the lights and music will stop.

Troubleshooting:

Problem: Mini UV Mirror Bead Chain does not play music or light up.

Action #1: Check the Batteries for strength and for proper battery polarity. Replace if weak or dead.

Problem: Mini UV Mirror Bead Chain does not work at all.

Action: #1 Check all Batteries for strength and for proper battery polarity. Replace if weak or dead.

Care of the Unit:

Mini UV Mirror Bead Chain can be wiped clean with any household multi-purpose, non-abrasive cleaner and disinfectant. We recommend Simple Green, which is a non toxic biodegradable all-purpose cleaner.

Do not use abrasive cleaners, as they will scratch the surface of the unit.

Do not submerge the unit, as it will damage the electrical components.