


Capability Switches Comparison Chart


Best Selling Switches

***Activation Force** - the amount of pressure needed on a switch to turn on a device. For your reference, a nickel weighs 5g.

DEVICE	PHOTO	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR
Gumball Switches #782R, Y, B, G, P, #7801		Auditory click, Tactile	25 g	Pressure	2½"D	3¼" x ¾" ¼ lb	Red, Yellow, Blue, Green, Purple
Glitter Gumball Switches #7830B, G, M, P, Y		Auditory click, Tactile	25 g	Pressure	2½"D	3¼" x ¾" ¼ lb	Red, Yellow, Blue, Green, Purple
Jumbo Switches #787R, Y, B, G, P, #789		Auditory click, Tactile	20 g	Pressure	5 1/8"	5 1/8" x 1" ¾ lb	Red, Yellow, Blue, Green, Purple
Mini-Gumball Switches #783R, Y, B, G, P, #7802		Auditory click, Tactile	35 g	Pressure	1¼"D	2" x ½" 1 oz	Red, Yellow, Blue, Green, Purple

Hand, Finger, Foot & Body Switches

DEVICE	PHOTO	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/WEIGHT	COLOR
Adjustable Angle Switch #1032		Auditory click, Tactile	25 g	Pressure	5"	6½" x 8" x 8" 1 lb	Yellow/Black
Button Click Switch #775		Auditory click, Tactile	40 g	Pressure	1 7/8"	2½" x 2½" ¼ lb	Yellow/Blue
Compact Switches (Set of 3) #745		Auditory click, Tactile	100 g	Pressure	¾" x 3/8"	¾" x 3/8" ¼ lb	Red, Blue, Yellow
Cushion Grip Switch (#3141 Small & #3137 Medium)		Tactile, Auditory Click	40 g	Pressure	5" x 2" x 2"	5" x 2" x 2" ¼ lb	Blue
Eye Blink Switch #9008		Visual	N/A	Eye Blink	N/A	4½" x 2½" x 1" ½ lb	Black/Silver
Finger Isolation Button #716		Tactile	100 g	Pressure	1½"	3" x 1½" 3 oz	Yellow

DEVICE	PHOTO	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR
Finger Switch #717		Tactile	15 g	Pressure		1½" x ½" 1 oz	Blue
Pinch Switch #800		Tactile	100 g	Pressure, Pinch	3" x ¼"	3" x ¼" 2 oz	Red
The Twitch Switch #1061		N/A	Adjustable from 1 g to 40 g	Pressure Adjustable	1½" x ½"	½ lb	White
Thumb Switch #802		Tactile	100 g	Pressure	1½"	4½" x 1" 3 oz	Yellow
Tilt Switch #990		None	N/A	Tilt Movement		3" x 2" x 1" 2 oz	Black
Wobble Switch #1070		Auditory click, Tactile	10 g	Pressure	4" Spring	24" x 2" ½ lb	Red/Yellow

Lights, Vibration & Music Switches

DEVICE	PHOTO	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR
Jumbo Switch with Vibration #791		Auditory click, Visual, Tactile. Vibration	20 g	Pressure	5"	5" x 1/2" 1 lb	Yellow
Jumbo Switch w/Lights, Music and Vibration #768		Auditory click, Visual, Tactile, Vibration Lights	20 g	Pressure, Visual	5"	5" x 1/2" 3/4 lb	Green
*Switch & Glow #894 *Not suitable for those who are allergic to Latex.		Tactile, Visual Light, Music Vibration, Lights	150 g	Pressure, Visual	N/A	7" x 6 1/2" x 7" 1 lb	Green/Blue
Textured Orbit Ball Switch #887		Auditory click, Visual, Tactile, Lights	90 g	Pressure	5 1/2" x 4 1/2"	5 1/2" x 4 1/2" 1 lb	Yellow/Blue

Switches for Visually Impaired							
DEVICE	PHOTO	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR
Bright Switch for the Visually Impaired #2045		Visual, Tactile Light	50 g	Pressure	N/A	7" x 3¾" 1 lb	Black
Bright Red Switch #262		Auditory click, Visual, Tactile, Lights & Vibration	50 g	Pressure	N/A	5" x 2½" ¾ lb	Blue
Tail Light Switch, Say It Play It #464		Lights. Visual, Tactile, Voice Output, Vibration, Auditory Click	90 g	Pressure	N/A	7¾" x 2½" 1 lb	Silver, Yellow

Pillow Switches							
DEVICE	PHOTO	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR
Pillow Switch on Easy Flex Mount #1014		Auditory click, Tactile	40 g	Pressure	3"	20" 1 lb	Black & Green
Petite Pillow Switches #831, S, R, G		Auditory click, Tactile	40 g	Pressure	3"	3" x 2" 3 oz	Green, Red, Safari, Floral
Plate Switches							
DEVICE	PHOTO	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR
Plate Switch #810		Tactile, Auditory Click	50 g	Pressure	3" x 5"	3" x 5" ½ lb	Yellow/Blue
Sensory Plate Switch #886		Auditory click, Visual, Tactile, Lights, Vibration	75 g	Pressure, Visual	5" x 8"	11¾" x 8½" x 3" 1¼ lbs	Yellow/Black

Saucer Switches							
DEVICE	PHOTO	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/WEIGHT	COLOR
Jumbo Switch with Latch Timer #792		Tactile	20 g	Pressure	5"	5" x 1/2" 1 lb	Red
Mini-Saucer Switch #718-S		Auditory click, Tactile	40 g	Pressure	3 1/4" x 3 1/4" x 1"	2 1/4" 1/4 lb	Blue/Yellow
Modified Pancake Switches #769R, B, Y, G, #766		Auditory click, Tactile	40 g	Pressure	5 3/8"	7" 1/2 lb	Red, Blue, Yellow, Green
On/Off Jumbo Switch #7805		Tactile	20 g	Pressure	5"	5" x 1/2" 1/2 lb	Yellow
Saucer Switch #718		Auditory click, Tactile	25 g	Pressure	7" x 1"	6" 1/2 lb	Blue/Yellow

Mouth, Sip, Puff & Breath Switches

DEVICE	PHOTO	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR
*Grip and Puff Switch #700		None	20 g	Pneumatic, Pressure or Puff	N/A	½ lb	White
*Sip and Puff Switch #970 *Not returnable due to sanitary reasons noted by OSHA standards. Straws should be sterilized before use.		None	20 g	Pneumatic, Pressure, Sip or Puff	Straw	19" Gooseneck 2¼ lbs	Black

Textured Switches

DEVICE	PHOTO	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR
Gooshy Switch w/Vibration, Lights & Music #636		Tactile, Visual, Lights, Music, Vibration, Auditory Click	90 g	Pressure	N/A	10" x 4½" x 2" ¼ lb	Blue/Yellow
*Koosh Switch #883 *Not suitable for those who are allergic to Latex.		Tactile, Visual Light, Music Vibration, Auditory Click	30 g	Pressure	N/A	8½" x 5¼" x 2½" 1 lb	Yellow
Multi-Sensory Center #725, #726		Auditory click, Tactile Light, Music Vibration	40 g	Pressure		9" x 6" 1 lb 11" x 7" 1¼ lbs	Blue/Yellow
Senso Dot Switch #682 & Set of 4 #683		Tactile, Vibration, Auditory Click	40 g	Pressure	6" x 2¾"	6" x 2¾" 1 lb	Yellow/Green, Color may vary
Small Oval Textured Switches (Blue) #711-S		Auditory click, Tactile	40 g	Auditory click, Tactile	4½"	5" x 4½" x 1" ¾ lb	Black/Yellow
Soft Shaggy Switch #3133		Tactile	40 g	Pressure	7"	7" x 5" x 2" ¼ lb	Green/Color may vary

DEVICE	PHOTO	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR
Super Bumpy Switch #736		Auditory click, Tactile	30 g	Auditory click, Tactile	8" x 4½"	8" x 4½" x 1/8" ¼ lb	Yellow
Textured Saucer Switches (Set of 3) #722		Auditory click, Tactile	25 g	Pressure	7"	6" 1¼ lbs	Orange/Blue, Blue/Yellow, Yellow/Red
Textured Roller Switch #639		Tactile, Visual, Lights, Vibration, Music	15 g	Pressure	N/A	11" x 4½" 1¼ lbs	Yellow/Green
Theme & Animal Switches							
DEVICE	PHOTO	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR
Lady Bug Switch w/Lights, Music & Vibration #729		Visual, Tactile, Lights, Vibration	50 g	Pressure	6" x 6¼" x 4"	9" x 6¼" x 4" 1 lb	Black/Red/ Yellow

Wheelchair & Bedside Switches

DEVICE	PHOTO	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR
Headswitch #705		Auditory click, Tactile	30 g	Pressure	3¼"	3¾" ½ lb	Black
Gumball Switch on Easy Flex Mount #1015		Auditory click, Tactile	30 g	Pressure	2½"	20" 1 lb	Black/Yellow
Pillow Switch on Easy Flex Mount #1014		Auditory click, Tactile	40 g	Pressure	3"	20" 1 lb	Black & Green
Ultimate Switch #1671		Auditory click, Tactile	10 g	Pressure	4" Spring	19" 2½ lbs	Black
Ultimate Switches #1551, 1552		Auditory click, Visual, Tactile	10 g	Pressure	4" Spring	6" 2 lbs 12" 2½ lbs	Black
Wobble Switch w/19" Gooseneck and Clamp #1063		Auditory click, Tactile	10 g	Pressure	4" Spring	19" 2¼ lbs	Black/Yellow

Wireless, Sound & Motion Activated Switches

DEVICE	PHOTO	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR
Movement Sensor Switch #655	 A small white rectangular device with a blue sensor arm and a black cable.	Visual	N/A	Motion, Adjustable sensitivity	N/A	4½" x 2½" x 1" ½ lb	White
Sound Activated Switch #1051	 A white rectangular device with a green indicator light and a black cable.		N/A	Sound	N/A	4½" x 2½" x 1" ¼ lb	White
Wireless Ultimate Switch #2700	 A person wearing a green shirt and a black head-mounted device with a sensor.	Auditory click, Tactile	10 g	Pressure	4" Spring	19" 2¾ lbs	Black

Switches for Computers and iPads

DEVICE	PHOTO	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR
Computer Mouse Interface #1165		Auditory click, Tactile, Visual	20.9	Switch Input Pressure	5"	5" ¾ lb	Yellow
iPad Wireless Switch #1164		Auditory click, Tactile, Visual	20.9	Switch Input Pressure	5"	5" ¾ lb	Yellow
iPad Wireless Two Switch #1167		Auditory click, Tactile, Visual	20.9	Switch Input Pressure	2 - 5"	5" x 12½" x ¾" ½ lb	Yellow/Red

Switch Modifiers

DEVICE	PHOTO	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR
Switch Modifier #605		N/A	N/A	Switch Input Pressure		4½" x 2½" x 1" ¼ lb	White

Switches Kits for Assessment							
DEVICE	PHOTO	FEEDBACK	ACTIVATION FORCE*	ACTIVATION TYPE	ACTIVATION SURFACE	SIZE/ WEIGHT	COLOR
Assessment Kit for More Involved Clients #7755		Varies	Varies	Pressure	Varies	Varies	Varies
Basic Switch Assessment Kit #7750		Varies	Varies	Pressure	Varies	Varies	Varies
Classroom Switch Kit #7803N		Auditory click, Tactile	Varies	Pressure	Varies	Varies	Varies
For further information please see the classroom kit section of our website.							